


How to Measure Bolster Indoor Cushion Cover

Measurement Instruction

- ✓ Give exact measurement from edge to edge.
- ✓ We add 0cm to -2.5cm size variance on the given width/depth for perfect fit.




Measurement (cm)

1. Diameter

2. Length


1. Diameter:

Diameter of the Bolster Cushion Cover



1. Length:

Length of the Bolster Cushion Cover



We encourage you to upload the image of the article you need the cover for - this helps our team ensure covers are designed keeping your requirements in mind.